

BATRIK

ACETABULAR REAMER BRUSH

Our double-ended cleaning brush is manufactured exclusively for cleaning acetabular bone reamers. The “Pig Tail” rounded style circular head brush is designed to effectively clean inside the reamer, while the rigid straight brush on the reverse side is used for cleaning individual cutting holes. Special fusion of the brush heads to the hand piece along with its high quality antimicrobial medical grade bristles ensures a durable, reusable, and highly effective cleaning tool.



STRAIGHT BRUSH



ROUND BRUSH



PRODUCT SELECTION

ITEM #	DESCRIPTION	HANDLE LENGTH	STRAIGHT BRUSH	ROUND BRUSH HEAD
50-ACETAB	Acetabular Reamer Cleaning Brush	4.72" / 120 mm Ø 0.47" / 12 mm	1.18" / 30 mm Ø 0.3125" / 8 mm	2" / 50 mm Ø 1.75" / 45 mm

PACKAGING: Bag of 3

NOTE: Specifications subject to small variances