

BATRIK

# AUTOCLAVE CLEANING BRUSH KITS



## ROUND BRUSH

- Features 3 holes for adjustable pole placement
- Adjustable brush placement (parallel or perpendicular) and longer brush length provides easy access to hard-to-reach corners
- Made with FDA approved polyester bristles

## HAND BRUSH

- Ergonomic handle for better hand comfort
- Longer, FDA approved polyester bristles
- Easy to clean (no cavity)

## PRODUCT SELECTION

### REGULAR AUTOCLAVE CLEANING BRUSH KIT

ITEM # (KIT)	DESCRIPTION	ITEM # (UNIT)	DESCRIPTION	DIMENSIONS
<b>50-CLAVEREG</b>	Autoclave REGULAR Brush Cleaning Kit	<b>50-35900</b>	1 Hand Brush	Length: 3.5" x 8.3" / 9 cm x 21 cm
		<b>50-55800</b>	1 Rounded Corner Brush (Head Brush)	Head size: 5.5 x 9.8" / 14 x 25 cm
		<b>50-EXTPOLE</b>	1 Anodized Aluminum Extension Pole	2-4 ft / 68-120 cm

### SUPERIOR AUTOCLAVE CLEANING BRUSH KIT

ITEM # (KIT)	DESCRIPTION	ITEM # (UNIT)	DESCRIPTION	DIMENSIONS
<b>50-CLAVESUP</b>	Autoclave SUPERIOR Brush Cleaning Kit	<b>50-35900</b>	1 Hand Brush	Length: 3.5" x 8.3" / 9 cm x 21 cm
		<b>50-55800</b>	1 Rounded Corner Brush (Head Brush)	Head size: 5.5 x 9.8" / 14 x 25 cm
		<b>50-EXTPOLE</b>	1 Anodized Aluminum Extension Pole	2-4 ft / 68-120 cm
		<b>SF-100</b>	2 SAVEFACE™ Face Shields	9" / 23 cm

## REPLACEMENT COMPONENTS



**HAND BRUSH**  
**50-35900**

Packaging: 1



**BRUSH HEAD**  
**50-55800**

Packaging: 1



**EXTENSION POLE**  
**50-EXTPOLE**

Packaging: 1



**SAVEFACE™**  
**FACE SHIELDS SF-100**

Packaging: 100 / Box