

BATRIK REUSABLE INSTRUMENT CLEANING BRUSHES


Designed specifically for hospital applications, our Reusable Instrument Cleaning Brushes feature long and flexible bristles to maximize the effectiveness of large-surface cleaning.

- High quality polyamide clear bristles
- Easy to identify debris
- Autoclavable
- Latex free

REUSABLE JAR/BOTTLE BRUSH



STERILIZER DRAIN BRUSH



Blue Tip
for added safety

Heavy Duty
Stainless Wire

Ergonomic
Handle

PRODUCT SELECTION

BLUE DELRIN PLASTIC

ITEM #	DESCRIPTION	HANDLE LENGTH	BRUSH LENGTH	BRUSH (Ø)	TOTAL LENGTH
50-10295	Reusable Jar / Bottle Brush	120 mm ± 2 mm	120 mm ± 5 mm	70 mm ± 2 mm	13.19" / 355 mm
50-10224	Sterilizer Drain Cleaning Brush	N/A	4" / 10.1 cm	2" / 5.08 cm	24" / 609 mm

PACKAGING: Bag of 3